

COHESIVE TAPE

LATEX or LATEX FREE COHESIVE TAPE 1/2"


PART CODE: 98062 (latex)

PART CODE: 98198 (latex free)

TAPE WIDTH: 1/2" ($\pm 1/32$ ") [~ 17.7mm (± 0.8)]

RECOMMENDED STORAGE: Store in cool dry place in original bag and box between bagged 60-80°F (15.5-26.7°C)

EXPIRATION DATE: 12 months from date of manufacture with proper storage



COHESIVE TAPE

LATEX or LATEX FREE COHESIVE TAPE 1"

PART CODE: 98216 (latex)

PART CODE: 98195 (latex free)

TAPE WIDTH: 1" ($\pm 1/32$ ") [$\sim 25.4\text{mm} (\pm 0.8)$]

RECOMMENDED STORAGE: Store in cool dry place in original bag and box between bagged 60-80°F (15.5-26.7°C)

EXPIRATION DATE: 12 months from date of manufacture with proper storage



ROHS. Restriction of the use of certain hazardous substances in electrical and electronic equipment which established a maximum concentration value up to 0.1% by weight for lead, mercury, chromium (VI), polybrominated biphenyls (PBB) or polybrominated diphenyl ethers (PBDE), bis (2-ethylhexyl)phthalate (DEHP), benzyl butyl phthalate (BBP), dibutyl phthalate (DIBP) and up to 0.01% for cadmium. Product DOES NOT contain any of the above listed substances

DESIGN AND SPECIFICATIONS ARE SUBJECT TO CHANGE WITHOUT PRIOR NOTICE. CODE 00635 REV. 00

www.tecnoidealsrl.com